FOR MORE AUTOMATION

UPGRADE NOW!






DentalCAD® 3.2 Elefsina at a glance

New features, more automation and seamless integration

- Preparation margin repair tool enables you to make the most out of imperfect intraoral scans: fix margins to get an acceptable result despite artifacts or holes in the scan
- Save time with Auto Articulator for automatic dynamic occlusion*
- Integration of compact milling machines:
 select your preferred block and start production
 seamlessly from within your CAD software*
- Switch tooth libraries while keeping previous setup






*Requires add-on module

Elevated Smile Creator for better visualization and pre-operative treatment planning

- Design smiles with virtual gingiva for superior outcomes
- Design smiles based on patients' natural teeth
- Improved 3D-to-photo alignment: adjust 3D view to match rendering to the camera's focal length



Bite Splint Module offers more design speed

- More **automation** in bite splint design
- Design bite splints for both arches simultaneously



Full Denture Module improved

- Copy Denture: Instead of relining an old denture, rapidly and inexpensively produce a new one
- Design dentures on ball attachments
- Add-on for the **Ivotion™ Denture System** now available*



*Ivotion library purchase from Ivoclar required

iTero-exocad Connector™ redefines seamless workflows

- Lab/doctor collaboration: doctors can now share external case-related files (images, videos and X-rays) via MyiTero™ portal
- New bi-directional case communication lets you share designs in 3D
- Validate margins with iTero intraoral camera images and visualization of tooth structure with iTero NIRI images



New beta feature: post and core*

- Design post and core restorations
- A crown matching the core can then be directly designed



*Available upon request for beta testers

Get the full highlights from your reseller.

exocad.com/dentalcad-elefsina

Find your exocad dealer here: exocad.com/partners



