












































# Approved Denture Tooth Libraries














987 mould chart libraries from 18 tooth manufacturers










Approved Denture Tooth Libraries					Design Options			Production Options					
	Included tooth setup presets	Included anterior upper moulds	Included anterior lower moulds	Included posterior moulds	 Scalable	 Bridged teeth	 Monoblock for Try-In	 Printed base & Printed teeth	 Printed base & Pre-manufactured teeth	 Printed base & Milled teeth	 Milled base & Pre-manufactured teeth	 Milled base & Trimmed pre-manufactured teeth	 Milled base & Milled teeth
 <b>CANDULOR AG</b> (56 moulds)													
<b>Bonartic® II NFC+</b>		0	0	3	●	●	●	●	●	●	●	●	●
<b>BonSelect TCR</b>		0	0	8	●	●	●	●	●	●	●	●	●
<b>Condyloform® II NFC+</b>		0	0	3	●	●	●	●	●	●	●	●	●
<b>PhysioSelect® TCR</b>		18	4	0	●	●	●	●	●	●	●	●	●
<b>PhysioStar® NFC+</b>		15	5	0	●	●	●	●	●	●	●	●	●
 <b>Dentsply Sirona Inc.</b> (115 moulds)													
<b>Genios®</b> including 1to1 and 1to2 posterior setups		17	10	18			●		●		●	●	(compatible with DS workflow)
<b>IPN 3D™ Portrait® Inspired</b> including lingualized and balanced setups		7	7	9			●		●		●	●	(compatible with DS workflow)
<b>Portrait® IPN®</b>		22	12	13			●		●		●	●	(compatible with DS workflow)














Approved Denture Tooth Libraries					Design Options			Production Options					
	Included tooth setup presets	Included anterior upper moulds	Included anterior lower moulds	Included posterior moulds	 Scalable	 Bridged teeth	 Monoblock for Try-In	 Printed base & Printed teeth	 Printed base & Pre-manufactured teeth	 Printed base & Milled teeth	 Milled base & Pre-manufactured teeth	 Milled base & Trimmed pre-manufactured teeth	 Milled base & Milled teeth
<b>exocad</b> exocad GmbH (18 moulds)													
generic denture	11	5	5	5	●	●	●	●	●	●	●	●	●
generic denture scalable automatic scaling	4	1	1	1	●	●	●	●	●	●	●	●	●
 gebdi GmbH (23 moulds)													
Tribos 501®		9	7	7	●	●	●	●	●	●	●	●	●
<b>HUGE</b> HUGE (164 moulds) (1/2)													
Bluebell		14	11	6	●	●	●	●	●	●	●	●	●
KAIJING		10	10	3	●	●	●	●	●	●	●	●	●
KAILI		9	5	4	●	●	●	●	●	●	●	●	●
KAIPLUS		9	9	3	●	●	●	●	●	●	●	●	●
MAIST		20	12	10	●	●	●	●	●	●	●	●	●

Approved Denture Tooth Libraries					Design Options			Production Options					
	Included tooth setup presets	Included anterior upper moulds	Included anterior lower moulds	Included posterior moulds	 Scalable	 Bridged teeth	 Monoblock for Try-In	 Printed base & Printed teeth	 Printed base & Pre-manufactured teeth	 Printed base & Milled teeth	 Milled base & Pre-manufactured teeth	 Milled base & Trimmed pre-manufactured teeth	 Milled base & Milled teeth
<b>HUGE</b> HUGE (2/2)													
Sonning		15	10	4	●	●	●	●	●	●	●	●	●
<b>ivoclar Ivoclar</b> (127 moulds)													
<b>Blueline®</b> activation needed: contact Ivoclar to unlock	56	28	7	6	●	●	●	●	●	●	●	●	●
<b>Ivotion®</b> activation needed: contact Ivoclar to unlock	36	19	3	6	●	●	●	●	●	●	●	●	●
<b>SR Phonares® II</b> activation needed: contact Ivoclar to unlock	52	18	6	5	●	●	●	●	●	●	●	●	●
<b>SR Vivodent® S DCL / S PE</b> activation needed: contact Ivoclar to unlock	48	16	7	6	●	●	●	●	●	●	●	●	●
<b>KULZER Kulzer GmbH</b> (99 moulds) (1/2)													
<b>Artic® digital</b>		21	9	5			●		●		●	●	
<b>Delara®</b>		18	5	4			●		●		●	●	
<b>Pala® Idealis</b>		0	0	4			●		●		●	●	

Approved Denture Tooth Libraries					Design Options			Production Options					
	Included tooth setup presets	Included anterior upper moulds	Included anterior lower moulds	Included posterior moulds	 Scalable	 Bridged teeth	 Monoblock for Try-In	 Printed base & Printed teeth	 Printed base & Pre-manufactured teeth	 Printed base & Milled teeth	 Milled base & Pre-manufactured teeth	 Milled base & Trimmed pre-manufactured teeth	 Milled base & Milled teeth
 <b>KULZER</b> <small>MITSUBI CHEMICALS GROUP</small> <b>Kulzer GmbH (2/2)</b>													
<b>Pala® Mondial</b>		18	10	5			●		●		●	●	
 <b>MAJOR</b> <b>MAJOR PRODOTTI DENTARI S.p.A. (23 moulds)</b>													
<b>MAJOR PLUS COMP</b>		10	4	9	●	●	●	●	●	●	●	●	●
 <b>Merz Dental GmbH (72 moulds)</b>													
<b>artegral® life Basic</b> contact Merz Dental for further info ● Free reduced data ○ Additional options with high resolution data (coming soon)	120	20	12	4	●	●	●	●	○ (available only with high resolution data)	●	○ (available only with high resolution data)	○ (available only with high resolution data)	●
 <b>myerson.</b> <b>Myerson LLC. (36 moulds)</b>													
<b>Kenson®</b>		18	9	9			●	●	●	●	●	●	●
Special (coming soon)													

Approved Denture Tooth Libraries					Design Options			Production Options					
	Included tooth setup presets	Included anterior upper moulds	Included anterior lower moulds	Included posterior moulds	 Scalable	 Bridged teeth	 Monoblock for Try-In	 Printed base & Printed teeth	 Printed base & Pre-manufactured teeth	 Printed base & Milled teeth	 Milled base & Pre-manufactured teeth	 Milled base & Trimmed pre-manufactured teeth	 Milled base & Milled teeth
 <b>Polident d.o.o.</b> (14 moulds)													
REF-LINE		4	1	9					●		●	●	
 <b>Ruthinium® - Dental Manufacturing S.p.A.</b> (80 moulds)													
Acry Plus including lingualized setups for all posterior moulds		28	12	10	●	●	●	●	●	●	●	●	●
Acry Smart		15	7	8	●	●	●	●	●	●	●	●	●
 <b>SHOFU® Inc.</b> (17 moulds)													
Veracia SA		9	4	4			●		●		●	●	
 <b>Unidesa Odi</b> (18 moulds)													
odident+		7	7	4	●	●	●	●	●	●	●	●	●

Approved Denture Tooth Libraries					Design Options			Production Options					
	Included tooth setup presets	Included anterior upper moulds	Included anterior lower moulds	Included posterior moulds	 Scalable	 Bridged teeth	 Monoblock for Try-In	 Printed base & Printed teeth	 Printed base & Pre-manufactured teeth	 Printed base & Milled teeth	 Milled base & Pre-manufactured teeth	 Milled base & Trimmed pre-manufactured teeth	 Milled base & Milled teeth
<b>VinciSmile® VinciSmile Group LLC</b> (40 moulds)													
<b>NOBILENT Prime DXL</b>		13	7	3	●	●	●	●	●	●	●	●	●
<b>NOBILENT Select XL</b>		7	6	4	●	●	●	●	●	●	●	●	●
<b>VITA VITA Zahnfabrik H. Rauter GmbH &amp; Co. KG</b> (28 moulds)													
<b>VIONIC VIGO®</b> including buccal, crossbite, lingualized and Dr. Gerber setups		8	4	16			●		●		●	●	
<b>YAMAHACHI DENTAL Yamahachi Dental Mfg.,Co.</b> (35 moulds)													
<b>CROWN PX &amp; NS</b>		8	6	0			●		●		●	●	
<b>EFUCERA A</b>		4	4	0			●		●		●	●	
<b>EFUCERA PX &amp; NS</b>		0	0	4			●		●		●	●	
<b>NAPERCE</b>		0	0	2			●		●		●	●	
<b>NEW ACE</b>		4	3	0			●		●		●	●	

Approved Denture Tooth Libraries					Design Options			Production Options					
	Included tooth setup presets	Included anterior upper moulds	Included anterior lower moulds	Included posterior moulds	 Scalable	 Bridged teeth	 Monoblock for Try-In	 Printed base & Printed teeth	 Printed base & Pre-manufactured teeth	 Printed base & Milled teeth	 Milled base & Pre-manufactured teeth	 Milled base & Trimmed pre-manufactured teeth	 Milled base & Milled teeth
 <b>HUXIAN YUWEI DENTAL MATERIAL CO.,LTD.</b> (22 moulds)													
Typlus		12	4	6	●	●	●	●	●	●	●	●	●
 <b>Davis Schottlander &amp; Davis Ltd.</b> (coming soon)													
enigmalive													
 <b>DEPRAG (Zhengzhou) Dental Technology Co., Ltd.</b> (coming soon)													
Euphon													
Euphon Pro													
Solidz MAX													
 <b>New Stetic S.A.</b> (coming soon)													
Bident													